

Table of contents

\MR-Physik

AG Zaiss

IMPI04

COMPARE_APTw_001_of_MPI04_816b_816c

[gre_cest_816B_APTw_001_FastTrec_3s](#)
[gre_cest_816C_APTw_001_slabsel_CS2](#)
[gre_cest_MPI04_APTw_001_FA8_fatsat](#)
[gre_cest_MPI04_WASAB1](#)

\\MR-Physik\AG Zaiss\IMPI04\COMPARE_APTw_001_of_MPI04_816b_816c\gre_cest_816B_APTw_001_FastTrec_3s

TA: 3:35 PM: FIX Voxel size: 2.0x2.0x5.0 mmPAT: 2 Rel. SNR: 1.00 : WIPCEST

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.1 ms
TE	2.00 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HC3-6

Contrast - Common

TR	4.1 ms
TE	2.00 ms
MTC	Off
Magn. preparation	None
Flip angle	8 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Multiple series	Off

Resolution - Common

FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
Slice oversampling	16.7 %

Geometry - Common

Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.1 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R1.2 A3.0 H19.2
R	1.2 mm
A	3.0 mm
H	19.2 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-6.5
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 A7.2 H20.4 mm
! Orientation	T > C-6.7
! Rotation	0.00 deg
! A >> P	220 mm

System - Adjust Volume

! R >> L	181 mm
! F >> H	80 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.257976 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	4.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	4.1 ms
Concatenations	1
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	220 mm
FoV phase	82.1 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	30
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s

Inline - Soft Tissue

Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8 deg
Measurements	30
Contrasts	1
TR	4.1 ms
TE	2.00 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Reordering	Spiral
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Bandwidth	700 Hz/Px

Sequence - Part 2

Define	Segments
Shots per slice	47
Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off

Sequence - Special

CEST	On
Pulse Type	Gauss
B1 definition	pure FA
B1	1.78 uT

Sequence - Special

Spiral Elong. E	0.5 (null)
No. of Pulses	36 n
Pulse Duration	49920 µs
Duration Scale Factor	1
Interpulse Delay	5 ms
Recover Time	3000 ms
Recover Time M0	3000 ms
Offset M0	300 ppm
Deactivate Duty Cycle Check	1
Spoiling	only last
Offset Distribution	File
Suppressed M0	On
Offset	4.0 ppm
Freq Shift	1.2 ppm
Delta Freq	1.2 ppm
Scale Factor	10.0 (null)
B0 Correction	PhaseMap
Image Registration	None
Save Unregistered	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\MR-Physik\AG Zaiss\IMPI04\COMPARE_APTw_001_of_MPI04_816b_816c\gre_cest_816C_APTw_001_slabsei_CS2

TA: 3:19 PM: FIX Voxel size: 2.0x2.0x5.0 mmPAT: Off Rel. SNR: 1.00 : fWIP

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.0 mm
TR	4.00 ms
TE	2.00 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HC3-6

Contrast - Common

TR	4.00 ms
TE	2.00 ms
Flip angle	8.0 deg
Fat suppr.	None
Lines Per Shot	496
Water suppr.	None
Dixon	Off
Dixon evaluation	Off

Contrast - Dynamic

Averages	1
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Multiple series	Off

Resolution - Common

FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.0 mm
Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
View sharing	Off
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm

Geometry - Common

FoV phase	82.1 %
Slice thickness	5.0 mm
TR	4.00 ms
Multi-slice mode	Sequential
Series	Ascending
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R1.2 A3.0 H19.2
R	1.2 mm
A	3.0 mm
H	19.2 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-6.5
> S	0.0

Geometry - Saturation

Fat suppr.	None
Water suppr.	None
Dixon	Off
Dixon evaluation	Off
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 A7.2 H20.4 mm
! Orientation	T > C-6.7
! Rotation	0.00 deg
! A >> P	220 mm
! R >> L	181 mm

System - Adjust Volume

! F >> H	80 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.257976 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.500
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

View sharing	Off
Flip angle	8.0 deg
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Burn time-to-center	On
Temporal interpolation	1
3D centric reordering	Off
Time to center	11.1 s

Inline - Inline

Subtract	Off
Measurements	30
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8.0 deg
Measurements	30
Contrasts	1
TR	4.00 ms
TE	2.00 ms

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	On
Asymmetric echo	Off
Contrasts	1
Optimization	None
Multi-slice mode	Sequential

Sequence - Part 1

Bandwidth	700 Hz/Px
-----------	-----------

Sequence - Part 2

RF pulse type	Normal
Gradient mode	Normal
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off

Sequence - Special

Sparse Sampling	On
Increment PE	1.0000
Increment Par	5.0000
Temporal Incoherence	0.0000
Reg. Param 1	0.0001
Reg. Param 2	0.0000
Additional scaling	0.0000
Acceleration Factor	2.0000
Channel Compression	1
ReconMode	2
Number of RefScansPE	36
Number of RefscansPar	12
Iterations	50
CSM Mode	2
Retro Recon Ini File	0
Max Power	5
Sparse Recon	Off
CEST	APT
Fill Time after FatSat Pulse	0 us

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\MR-Physik\AG Zaiss\IMPIO4\COMPARE_APTw_001_of_MPI04_816b_816c\gre_cest_MPI04_APTw_001_FA8_fatsat

TA: 3:16 PM: FIX Voxel size: 2.0×2.0×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl_CEST

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.0 ms
TE	2.00 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HC3-6

Contrast - Common

TR	4.0 ms
TE	2.00 ms
MTC	Off
Magn. preparation	None
Flip angle	8-deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s
Multiple series	Off

Resolution - Common

FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
Slice oversampling	16.7 %

Geometry - Common

Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.0 ms
Multi-slice mode	Interleaved
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R1.2 A3.0 H19.2
R	1.2 mm
A	3.0 mm
H	19.2 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-6.5
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 A7.2 H20.4 mm
! Orientation	T > C-6.7
! Rotation	0.00 deg
! A >> P	220 mm

System - Adjust Volume

! R >> L	181 mm
! F >> H	80 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.257976 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	4.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	4.0 ms
Concatenations	1
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	220 mm
FoV phase	82.1 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	30
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	30
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s

Inline - Soft Tissue

Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Pause after meas. 24	0.0 s
Pause after meas. 25	0.0 s
Pause after meas. 26	0.0 s
Pause after meas. 27	0.0 s
Pause after meas. 28	0.0 s
Pause after meas. 29	0.0 s

Sequence - Special

No. of Pulses	36 n
Pulse Duration	49920 µs
Interpulse Delay	5 ms
Recover Time	2400 ms
Recover Time M0	2400 ms
No. of Locking Pulses	1 n
Lock Interpulse Delay	1 ms
Spoiling	only last
Fat Sat	On
Spiral Elong. E	0.5
Offset Distribution	APT
Suppressed M0	On
Offset	3.0 ppm
Freq Shift	0.0 ppm
Delta Freq	0.0 ppm
Scale Factor	10.0 (null)
Duration	1.0 ms
BWTP	9.6
EmpFactor	0.87
Samples	200
B0 Correction	Z-Spectrum
Image Registration	None

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Assistant

Mode	Off
Allowed delay	0 s

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	8-deg
Measurements	30
Contrasts	1
TR	4.0 ms
TE	2.00 ms

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Reordering	Spiral
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Bandwidth	700 Hz/Px

Sequence - Part 2

Define	Segments
Shots per slice	47
Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off

Sequence - Special

CEST	On
Pulse Type	Gauss
B1 definition	pure FA
B1	1.78 uT

\\MR-Physik\AG Zaiss\IMPI04\COMPARE_APTw_001_of_MPI04_816b_816c\gre_cest_MPI04_WASAB
1
TA: 2:02 PM: FIX Voxel size: 2.0×2.0×5.0 mmPAT: 2 Rel. SNR: 1.00 : fl_CEST

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Phase oversampling	0 %
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.0 ms
TE	2.00 ms
Averages	1
Concatenations	1
Filter	Raw filter
Coil elements	HC3-6

Contrast - Common

TR	4.0 ms
TE	2.00 ms
MTC	Off
Magn. preparation	None
Flip angle	6 deg
Fat suppr.	None
Water suppr.	None
SWI	Off

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	24
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s

Contrast - Dynamic

Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s
Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s
Multiple series	Off

Resolution - Common

FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
Base resolution	112
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off

Resolution - iPAT

PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	32
Accel. factor 3D	1
Ref. lines 3D	12
Reference scan mode	GRE/separate

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	On
Elliptical filter	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
Slice oversampling	16.7 %
Slices per slab	12
FoV read	220 mm
FoV phase	82.1 %
Slice thickness	5.00 mm
TR	4.0 ms
Multi-slice mode	Interleaved

Geometry - Common

Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	R1.2 A3.0 H19.2 mm
Orientation	T > C-6.5
Phase enc. dir.	R >> L
AutoAlign	---
Initial Position	R1.2 A3.0 H19.2
R	1.2 mm
A	3.0 mm
H	19.2 mm
Initial Rotation	90.00 deg
Initial Orientation	T > C
T > C	-6.5
> S	0.0

Geometry - Saturation

Saturation mode	Standard
Fat suppr.	None
Water suppr.	None
Special sat.	None

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	0 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
AutoAlign	---
Coil Select Mode	On - AutoCoilSelect

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L0.0 A7.2 H20.4 mm
! Orientation	T > C-6.7
! Rotation	0.00 deg
! A >> P	220 mm
! R >> L	181 mm
! F >> H	80 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
Excitation	Slab-sel.

System - Tx/Rx

Frequency 1H	123.257976 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	4.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	4.0 ms
Concatenations	1
Segments	1

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off
FoV read	220 mm
FoV phase	82.1 %
Phase resolution	100 %

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	24
StdDev	Off
Liver registration	Off
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Soft Tissue

Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off
Measurements	24
Pause after meas. 1	0.0 s
Pause after meas. 2	0.0 s
Pause after meas. 3	0.0 s
Pause after meas. 4	0.0 s
Pause after meas. 5	0.0 s
Pause after meas. 6	0.0 s
Pause after meas. 7	0.0 s
Pause after meas. 8	0.0 s
Pause after meas. 9	0.0 s
Pause after meas. 10	0.0 s
Pause after meas. 11	0.0 s
Pause after meas. 12	0.0 s

Inline - Soft Tissue

Pause after meas. 13	0.0 s
Pause after meas. 14	0.0 s
Pause after meas. 15	0.0 s
Pause after meas. 16	0.0 s
Pause after meas. 17	0.0 s
Pause after meas. 18	0.0 s
Pause after meas. 19	0.0 s
Pause after meas. 20	0.0 s
Pause after meas. 21	0.0 s
Pause after meas. 22	0.0 s
Pause after meas. 23	0.0 s

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Inline - MapIt

Save original images	On
MapIt	None
Flip angle	6 deg
Measurements	24
Contrasts	1
TR	4.0 ms
TE	2.00 ms

Sequence - Special

Suppressed M0	On
Offset	1.8 ppm
Freq Shift	0.0 ppm
Delta Freq	0.0 ppm
Scale Factor	10.0 (null)
Duration	1.0 ms
BWTP	9.6
EmpFactor	0.87
Samples	200
B0 Correction	Z-Spectrum
Image Registration	None

Sequence - Assistant

Mode	Off
Allowed delay	0 s

Sequence - Part 1

Introduction	On
Dimension	3D
Elliptical scanning	On
Phase stabilisation	Off
Reordering	Spiral
Asymmetric echo	Off
Contrasts	1
Flow comp.	No
Multi-slice mode	Interleaved
Bandwidth	700 Hz/Px

Sequence - Part 2

Define	Segments
Shots per slice	47
Segments	1
Acoustic noise reduction	None
RF pulse type	Fast
Gradient mode	Fast
Excitation	Slab-sel.
RF spoiling	On
Incr. Gradient spoiling	Off

Sequence - Special

CEST	On
Pulse Type	Rect
B1 definition	pure FA
B1	3.70 uT
No. of Pulses	1 n
Pulse Duration	5120 µs
Interpulse Delay	100 ms
Recover Time	2500 ms
Recover Time M0	12000 ms
No. of Locking Pulses	1 n
Lock Interpulse Delay	1 ms
Spoiling	only last
Fat Sat	On
Spiral Elong. E	0.5
Offset Distribution	Regular